



## TECHNICAL DATA

Definition/Water quality	Drinking water (water for human consumption)	Feed water for reverse osmosis	Dialysis water/permeate (water for diluting concentrated Haemodialysis solutions)		
		GLT			
Directive	98/83/EC	98/83/EC + procedural limit values	ISO 13959	European Pharmacopoeia	Recommendation applied to hygiene.
Chemical/physical parameters [ppm]					
Sodium (Na)	200	200	70	50	50
Potassium (K)		--	8	2	2
Calcium (Ca)		Total hardness	2	2	2
Magnesium (Mg)		< 1°dH or < 1.79°f	4	2	4
Boron (B)	1.0	1			
Barium (Ba)		0.7	0.1		0.1
Beryllium (Be)		0.004	0.0004		0.0004
Ammonium (NH <sub>4</sub> )	0.5	0.1		0.2	0.2
Aluminium (Al)	0.2	< 0.01	0.01	0.01	0.01
Metals					
– Copper (Cu)	2	1	0.1	--	0.1
– Arsenic (As)	0.01	0.01	0.005	--	0.005
– Lead (Pb)	0.01	0.01	0.005	--	0.005
– Silver (Ag)	-- 0.05	0.1	0.005	--	0.005
– Chromium (Cr)	0.01	0.05	0.014	--	0.014
– Selenium (Se)	0.005	0.01	0.09	--	0.01
– Stibium (Sb)	0.001	0.005	0.006	-- 0.001	0.005
– Mercury (Hg)	0.02	0.001	0.0002	--	0.0002
– Nickel (Ni)	-- 0.2	0.02	--	--	--
– Tin (Sn)	0.005	--	--	--	--
– Iron (Fe)	-- 0.05	< 0.1	-- 0.001	-- 0.1	-- 0.001
– Cadmium (Cd)	0.010	0.005	0.1	--	0.1
– Zinc (Zn)	--	5.0	--	--	--
– Manganese (Mn)		< 0.01	-- 0.002	-- 0.1	--
– Uranium (U)		0.01	0.1		--
– Thallium (Ti)		--			
or the sum of heavy metals					
Cyanide (CN)	0.05	0.05			0.02
Chlorine (Cl <sub>2</sub> )		Total chlorine: 0.0	0.1	0.1	0.1
1.2-dichlorethane	0.0030				
Chloramine					0.1
Chloride (Cl)	250	250		50	50
Fluoride (F)	1. May	1. May	0.2	0.2	0.2
Sulphate (SO <sub>4</sub> )	250	240	100	50	50

Definition/Water quality	Drinking water (water for humans consumption)	Feed water for Reverse Osmosis <i>GLT HD RO REVERSE OSMOSIS</i>	Dialysis water/permeate (water for diluting concentrated Hemodialysis solutions)		
			ISO 13959	European Pharmacopoeia	Recommendation applied hygiene
Directive	98/83/EC	98/83/EC + procedural limit values	ISO 13959	European Pharmacopoeia	Recommendation applied hygiene
Nitrate (NO <sub>3</sub> )	50	10	2 (as N)	2	2
Nitrite (NO <sub>2</sub> )	0.5	0.5			
Polycyclic aromatic hydrocarbons	0.00010	0.0001			
Benzene	0.0010	<b>0.001</b>			
Bromate	0.010	<b>0.01</b>			
Tetrachlorethene and trichlorethene	0.010	<b>0.005</b>			
Trihalogenmethane	0.050	<b>0.05</b>			
Vinyl chloride	0.00050	<b>0.0005</b>			
Silicic acid (SiO <sub>2</sub> )		<b>&lt; 10</b>			
pH value	6.5 – 9.5	<b>6.5 – 9.0</b>			
Temperature		<b>6 – 30 °C</b>			
Spec. conductivity	2500 µS/cm at 20 °C	<b>&lt; 1000 µS/cm at 20 °C</b>			
Silting index SDI <sub>(15)</sub> Clouding (NTU)	NTU < 1	<b>SDI (15 min) &lt; 5 (GLT RO) &lt; 3 (GLT RO) As per ASTM 4189</b>			
<b>Microbiological parameters</b>					
Total germ count] [CFU/ml]	< 100 (22 ± 2 °C, 44 ± 4h) < 100 (36 ± 1 °C, 44 ± 4h)	< 100 (22 °C) < 100 (36 °C)	< 100 (action at 50%) (17–23 °C, 7d)	< 10 <sup>2</sup> (30–35 °C, 5 d)	< 100 acc.to RKI (22 ± 2 °C, 3–7 d)
Enterococci	0 CFU/100ml	0 CFU/100ml			
E.-Coli/ coliform	0 CFU/100ml	0 CFU/100ml			
Endotoxins [EU/ml]			<0, 25 (action at 50%)	< 0.25	<0.25

a "Guideline for applied hygiene in dialysis units", ISBN 978-3-00-044348-0, 2013

**Note:**

Directive 98/83/EC and ISO 13959 specify limit values for rare substances that are not listed here; these can be looked up in the original publications. Compared to earlier publications, no information regarding phosphate is provided.

## Specifications

<b>Hemodialysis Device</b>	Up to 50 machines running at 800 mL/h per 2000 L/h
<b>Permeate Capacity</b>	Up to 75%
<b>Efficiency/Yield Dimensions</b>	180 x 95 x 258 cm
<b>(h x w x d) Weight (filled)</b>	400 kg
<b>Concentrate Pressure</b>	Max. 25 BAR

## Electrical Supply

<b>Electrical Supply/Three-phase</b>	380 V 3/N/PE, 50 Hz, 8 kW
<b>Overcurrent Protection</b>	32 A tripping characteristic (depending on voltage/version) D or K or similar recommended (due to high starting currents)
<b>Socket</b>	380 V: hardwired
<b>Type of Protection Against Electric Shock</b>	Protection Class I Type
<b>Applied Parts Classification</b>	B
<b>Degree of Ingress Protection Against Liquids</b>	Drip-proof
<b>Leakage Currents</b>	According to EN 60601-1 II
<b>Overvoltage Category</b>	II
<b>Material Group</b>	III b
<b>Operating Mode</b>	Continuous operation (standby)

## Product Water Quality

<b>Bacteria (CFU) and Endotoxins (EU)</b>	>99%
Product water quality depends on inlet water quality	>96%

<b>Feed Pressure</b>	Dynamic 2–6 BAR
<b>Minimum Inlet</b> Minimum inlet flow in liters per hour at maximum outlet capacity and a yield of 75%	<b>GLT HD 50 Plus + DP: min. 2500 L/h</b>
<b>Permeate Connection Inlet</b>	
<b>Water Connection</b>	Direct PE-Xa connector 25x3.5 (feed and return) on the system 1" external thread
<b>Drain Water Connection</b>	DN 32 (HT pipe)

## Operating Conditions

<b>Water Hardness</b>	<1.0 °dH
<b>Iron</b>	<0.1
<b>Manganese</b>	<0.1
<b>Chloride Silicate</b>	<100
<b>Total Chlorine</b>	<25 mg/L
<b>Feed Water Conductivity</b>	0.1 mg/L
<b>Total Salt Content</b>	<2500 uS/cm
<b>pH</b>	1500 mg/L
<b>Silt Density</b>	6–8
<b>Feed Water Temperature</b>	< Min. 5°C/max. 35°C Ambient pressure: 700–1150 hPa
<b>Ambient Temperature Range</b>	+5°C to +35°C
<b>Relative Humidity</b>	Up to 80% at 20°C (non-condensing)

## Transport and Storage Conditions

<b>Storage Temperature Range</b>	+5°C to +40°C (protect from freezing)
<b>Storage Time</b>	Storage time of preserved system: maximum 12
<b>Atmospheric Pressure</b>	Ambient pressure: 500–1150 hPa
<b>Relative Humidity</b>	Up to 80% at 20°C (non-condensing)

## Filling volumes of preservative / antifreeze

Number of modules	Sodium metabisulfite [gr]	MgCl [gr]	Glycerin 86% for -5°C [litres]	Glycerin 86% for -9°C [litres]	Glycerin 86% for -17°C [litres]	Total volume of liquid for RO [litres]
1	450	30	9.0	12.8	18.0	90
2	550	35	11.0	15.7	22.0	110
3	650	40	13.0	18.5	26.0	130
4	750	45	15.0	21.5	30.0	150
5	850	50	17.0	24.5	34.0	170
6	950	55	19.0	27.5	38.0	190

### Preservation using sodium metabisulfite

- If biofouling of the membrane is not to be expected and if the membrane is to be protected for storage, a solution with 0.5 wt./vol.% Sodium metabisulfite can be used.  
9.5% wt./vol.% glycerin must be added to guarantee frost protection down to -5 °C.
- It is helpful to produce a basic solution with sodium metabisulfite in a 20-fold concentration and to fill the supply tank with this basic solution.

### Stabilization

- 200 – 350 mg/l magnesium chloride (in the form of MgCl<sub>2</sub>) must be added to maintain membrane stability if this solution is to be stored for longer than one month.
- Let the preservative solution circulate through the membrane. Recirculate the solution through the mixing tank for one hour. The temperature must not exceed 35 °C.
- It is helpful here, too, to produce a basic solution with MgCl<sub>2</sub> in a 20-fold concentration and to fill the supply tank with this basic solution.

## OPERATING INDICATIONS

GLT systems consistently deliver **ultrapure water**, minimizing risks from sodium, potassium, calcium, aluminium, heavy metals, chlorine, and endotoxins. This makes GLT the trusted solution for patient safety.

### Note:

(\*\*optional) These features are not included in the standard version and are available only upon demand.  
(customized by project)



### Main Address

Ehlibeyt Mh., Tekstilciler Cd., No: 35/7,  
Çankaya, Ankara

### Istanbul Address

Girne Mh., Güneşlikli Sk., No: 9/A,  
Maltepe, Istanbul  
Fax: +90 312 247 50 25

Email: [info@galatamedical.com](mailto:info@galatamedical.com)  
[www.galatamedical.com](http://www.galatamedical.com)

When ordering, please include the following information:

- VAT Number Hospital or Dialysis Center Name
- Complete Delivery Address Requested Delivery
- Date Purchase Order Number Purchase
- New or Refurbished
- Contact Name and Contact
- Telephone Number E-mail address



EN ISO 13485:2016

