

GLT TEMPORARY HEMODIALYSIS CATHETER KIT



CATHETER

- ✓ Dual-D lumen
- ✓ Single dilator with gradually varied O.D.
- ✓ Rotatable suture hub
- ✓ Flexible polyurethane material
- ✓ Latex-free

GLT TEMPORARY HEMODIALYSIS CATHETER KIT



Standard Package:

Hemodialysis catheter (1)
 Guidewire (1)
 Introducer Needle (0/1)
 Y-Introducer needle (0/1)
 Dilator (1)
 Blue Introducer Syringe (0/1)
 Heparin Cap (2/3)
 Syringe (1/2)
 Needle (1)
 Scalpel (1)

Catalogue No.	Lumen Type	Catheter Type	Catheter O.D.	Catheter Length (cm)
GLTBSJ(Y)2-115-13(D)	Double-Lumen	Straight	11.5F	13cm
GLTBCJ(Y)2-115-13(D)	Double-Lumen	Curved	11.5F	13cm
GLTBSJ(Y)2-115-16(D)	Double-Lumen	Straight	11.5F	16cm
GLTBCJ(Y)2-115-16(D)	Double-Lumen	Curved	11.5F	16cm
GLTBSJ(Y)2-115-20(D)	Double-Lumen	Straight	11.5F	20cm
GLTBCJ(Y)2-115-20(D)	Double-Lumen	Curved	11.5F	20cm
GLTBSJ(Y)2-12-13(D)	Double-Lumen	Straight	12F	13cm
GLTBCJ(Y)2-12-13(D)	Double-Lumen	Curved	12F	13cm
GLTBSJ(Y)2-12-16(D)	Double-Lumen	Straight	12F	16cm
GLTBCJ(Y)2-12-16(D)	Double-Lumen	Curved	12F	16cm
GLTBSJ(Y)2-12-20(D)	Double-Lumen	Straight	12F	20cm
GLTBCJ(Y)2-12-20(D)	Double-Lumen	Curved	12F	20cm
GLTBSJ(Y)3-12-13(D)	Triple-Lumen	Straight	12F	13cm
GLTBCJ(Y)3-12-13(D)	Triple-Lumen	Curved	12F	13cm
GLTBSJ(Y)3-12-16(D)	Triple-Lumen	Straight	12F	16cm
GLTBCJ(Y)3-12-16(D)	Triple-Lumen	Curved	12F	16cm
GLTBSJ(Y)3-12-20(D)	Triple-Lumen	Straight	12F	20cm
GLTBCJ(Y)3-12-20(D)	Triple-Lumen	Curved	12F	20cm

