


GLT AV BLOODLINE TUBING SET



Optimized for high-volume dialysis centers



- ✓ **Sterile, single-use**
- ✓ **Medical-grade PVC**
- ✓ **Stable blood flow**
- ✓ **Flexible tubing**
- ✓ **Leak-proof connections**
- ✓ **Easy handling**

WHY GLT?

- Consistent product quality
- Competitive pricing
- Flexible OEM solutions
- Fast & stable supply

COMPATIBILITY & OPTIONS

- Compatible with major dialysis systems
- Custom configurations available
- CE-compliant production
- Available in DEHP-free option
- Designed for consistent daily use

Reliable Quality. Competitive Pricing. Trusted Supply.

info@galatamedical.com | www.galatamedical.com

CE-marked medical device

GLT AV BLOODLINE SETS

Designed for reliable blood circulation during dialysis treatment

TECHNICAL SPECIFICATIONS

Sterilization: ETO sterilized
Material: Medical-grade PVC
Tube Structure: Arterial & venous lines
Pump Seg.: Compatible with peristaltic pump
Connection Type: Luer lock
Drip Chamber: With/without filter
Injection Port: Available
Designed for patient safety
Disposable: Single-use
Pyrogen-free
Shelf life: 3 years

Model	Compatibility	Air Chamber
20 mL	Standard dialysis systems	20 mL
30 mL	Fresenius, Gambro, Nipro, B.Braun	30 mL
50 mL	Fresenius (large), Gambro	50 mL

Technical Advantages

Smooth inner surface
Stable performance
Reduced flow resistance design
Kink-resistant tubing structure

Standard Bloodline Sets

Wide range of model configurations
Compatible with major dialysis systems
Available in multiple chamber volumes
Clamps Color-coded (Red/Blue)
With Drain Bag (1.5L- 2.0L)
Standard (Included in Pack)

MODELS & CONFIGURATIONS

Model Range

HDI / HDJ / HDK / HDL / HDM / HDN / HDB Series

Available volumes:

20 mL 30 mL 50 mL 60 mL 25 mL

DEHP-Free Options

HDTA / HDTF / HDTZ Series

DEHP-free options are usually available upon request
Available for sensitive patient groups
Same performance with safer material option
Ensuring maximum patient safety

HDF ACCESSORY TUBING

Type	Ref & Model	Pump Tube Length	Compatibility
Standard	GLT-HDTUBE-030-067W	6.4&9.8mm 146x	Fresenius, B.Braun, Gambro, Nipro, Nikkiso
Unique	HDIT-02	430 mm	Fresenius 5008

Reliable Quality. Competitive Pricing. Trusted Supply.
info@galatamedical.com | www.galatamedical.com



CE-marked medical device