



HEMODIALYSIS CONSUMABLES

Water treatment system is the heart of hemodialysis treatment process. Our aim is to protect the heart of your system with our EU standart designs and solutions.

## REVERSE OSMOSIS SYSTEM

- ✓ Adjustable alarm limits
- ✓ Made in Türkiye
- ✓ Special type RO membranes
- ✓ Highest efficiency
- ✓ High pressure pumps made of 316 L stainless steel



## Standart Equipments and Specifications

- Stainless steel frame
- Conductivity meter
- Temperature meter
- Permeate and reject water flowmeters
- Water inlet, RO membrane inlet and outlet manometers
- Low pressure switches



## Technical Specifications

Model	Capacity (1 / day)	Sizes (LxWxH) (mm)
HD 1 PLUS	1 HD	750mm x 450mm x 850mm
HD 2 PLUS	2 HD	750mm x 450mm x 850mm
HD 5 PLUS	5 HD	640mm x 750mm x 1650mm
HD 10 PLUS	10 HD	640mm x 750mm x 1650mm
HD 15 PLUS	15 HD	640mm x 750mm x 1650mm
HD 20 PLUS	20 HD	640mm x 750mm x 1650mm
HD 30 PLUS	30 HD	800mm x 750mm x 1650mm
HD 40 PLUS	40 HD	800mm x 2400mm x 1800mm
HD 50 PLUS	50 HD	1000mm x 2400mm x 1800mm
HD 60 PLUS	60 HD	1000mm x 2400mm x 1800mm
HD 70 PLUS	70 HD	1000mm x 2400mm x 1800mm
HD 80 PLUS	80 HD	1000mm x 2400mm x 1800mm
HD 90 PLUS	90 HD	1000mm x 2400mm x 1800mm
HD 100 PLUS	100 HD	1000mm x 2400mm x 1800mm

